4R104=60, 4R104=75

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	26
PT	Manual de utilização	32
NL	Gebruiksaanwijzing	38
SV	Bruksanvisning	44
DA	Brugsanvisning	49
NO	Bruksanvisning	55
FI	Käyttöohje	60
PL	Instrukcja użytkowania	66
HU	Használati utasítás	72
CS	Návod k použití	77
RO	Instrucțiuni de utilizare	83
HR	Upute za uporabu	88
SL	Navodila za uporabo	94
SK	Návod na používanie	99
BG	Инструкция за употреба	105
TR	Kullanma talimatı	111
EL	Οδηγίες χρήσης	116
RU	Руководство по применению	123
JA	取扱説明書	129
ZH	使用说明书	135
KO	사용 설명서	139


1


2

3

4

INFORMATION

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INFORMATION

Last update: 2016-05-02

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1.1 Construction and Function

The 4R104=60 and 4R104=75 double adapters, sliding, are used to connect 2 pyramid adapters. The pyramid receivers make it possible to make distal and proximal angle adjustments of the prosthetic components in the frontal and sagittal planes. The dovetail guide makes it possible to shift the prosthetic components in the frontal and sagittal planes.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other components of this modular system. Any restrictions regarding the combination possibilities are listed in this section. Contact the manufacturer with questions.

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Approved for a body weight of **up to 100 kg**.

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the patient's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
---	---

3.2 General safety instructions

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified field of application (see Page 9).

CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service life").
- ▶ Only use the product for a single patient.

CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.

- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Double adapter, sliding	–
	Pre-assembled:	
8	Set screw	506G3=M8x16
1	Set screw	506G3=M5x8
2	Cap screw	501Z16

5 Preparation for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Installation in the modular prosthesis

Adjusting the stop

The set screw in the dovetail guide limits sliding of the adapter in one direction. For sliding in the opposite direction, the adapter can be turned or the set screw can be moved.

> **Required tools:** 710D1 torque wrench

► **CAUTION! Only use the adapter with the set screw screwed into the dovetail guide.**

Make sure that the set screw (see fig. 4) is screwed correctly into the dovetail guide (tightening torque: **5 Nm**).

Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

> **Required materials:** 710D4 torque wrench, 636K13 Loctite® 241

1) **Trial fitting:**

Screw in the set screws.

Use the torque wrench to tighten the set screws (**10 Nm**).

2) **Definitive mounting:**

Use Loctite® to secure the set screws.

Screw in the set screws.

Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).

3) Replace any set screws that are protruding or recessed too much with suitable ones (see selection table).

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

5.2 Sliding the adapter

> **Required tools:** 710D4 torque wrench

> **Prerequisite:** No load on the prosthesis and the prosthesis is not being worn.

1) Loosen the cap screws of the dovetail guide (see fig. 2).

2) Slide the adapter to the desired position.

- 3) Tighten the cap screws (installation torque: **15 Nm**).

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH