




Instructions for Use







D/P FLEXION™

Partnumber: K2D0xyyz(BR), FBD0xyyz (BR)

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	16
NO	Bruksanvisning	19
DA	Brugsanvisning	21
SV	Bruksanvisning	23
FI	Käyttöohjeet	25
NL	Gebruiksaanwijzing	28
PT	Instruções de Utilização	31
JA	取扱説明書	34
KO	사용 설명서	36

ENGLISH

This manual is intended for use by a certified Prosthetist.

D/P FLEXION™

- Maximum user weight 136 kg (299 lbs)*
- Only to be used with K2 Sensation or Flex-Foot Balance*
- Only for low activity levels (K2)*
- * Warranty voided if this requirement is not respected.

MEDICAL INDICATION

Lower extremity amputations, for example due to:

- Trauma
- Vascular Disease
- Cancer
- Congenital Defects

CONTRAINDICATIONS FOR USE

The foot should not be used in any sport activities.

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

Weight Kg	36-44	45-61	62-80	81-97	98-114	115-136
Weight lbs	79-97	98-134	135-176	177-214	215-251	252-300
21	0	1	N/A	N/A	N/A	N/A
22	0	1	N/A	N/A	N/A	N/A
23	0	1	2	N/A	N/A	N/A
24	0	1	2	N/A	N/A	N/A
25	N/A	1	2	3	4	5
26	N/A	1	2	3	4	5
27	N/A	N/A	2	3	4	5
28	N/A	N/A	2	3	4	5

INSTRUCTIONS FOR USE

The foot is a non-invasive, single patient, reusable prosthetic device.

- Bench- and dynamic alignment with alignment inserts. **(Figure 3a)**
- Inserts and O-rings shall be removed after alignment. **(Figure 3b)**
- After alignment take out alignment inserts. **(Figure 4)**
- If there is no movement after removing inserts while walking, lower the setting for both the D-valve and P-valve. **(Figure 5)**
- Start adjusting the P-valve **(Figure 6)**
 - If the heel is too stiff, decrease the P-valve setting.
 - If the heel is too soft, increase the P-valve setting.
- Start adjusting the D-valve **(Figure 7)**
 - If motion in stance is too slow, decrease the D-valve setting.
 - If motion in stance is too fast, increase the D-valve setting.
- Factory settings are set to 3 for both valves. **(Figure 8)**

CARE & CLEANING

Do not expose adapter to extreme temperature, sand, water or solvents. Should any exposure occur, clean adapter with a soft, damp cloth. Warranty is voided if these requirements are not respected.

RETROFITTING OF D/P FLEXION TO K2 SENSATION OR FLEX-FOOT BALANCE

In cases where D/P Flexion is being retrofitted to K2 Sensation or Flex-Foot Balance, or any other case where the D/P Flexion module needs to be attached to the foot the following torque specifications apply: 15Nm. Loctite 243 also needs to be applied to threads.

CAUTION

Avoid placing hands or fingers near moving parts. Please be aware that D/P Flexion is nickel coated and is filled with hydraulic fluid. Both of which can cause allergic reaction if it comes in contact with bare skin.

USE ENVIRONMENT

- Environmental operating temperature: 0°C (+32°F) -40°C (+105°F).
- Operating humidity: 10%-95% Relative humidity.
- Shipping and storage temperature: -15°C (+5°F) to 50°C (+122°F).

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CE CONFORMITY

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guideline.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee’s activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks. In the standard mentioned, test

levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads have been adapted adequately on the basis of the specified load level.

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceed!

For specific conditions and limitations of use see manufacturer’s written instructions on intended use!

ISO 10328 - label		
Category	Weight (Kg)	Label text
1	24	ISO 10328–P3–24–kg
2	34	ISO 10328–P3–34–kg
3	45	ISO 10328–P3–45–kg
4	55	ISO 10328–P3–55–kg

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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