

1. Open knee immobilizer and remove stay assemblies.
2. Position knee immobilizer so patella strap will go over the mid-line of the patella with wide end toward thigh and back stays along midline of leg.
3. Use small hook tab to hold body around knee.
4. Attach stay assemblies, inside first, to accommodate calf and thigh circumference.
5. Starting at the bottom and moving up, thread straps through buckles and pull back, attaching straps where comfortable.
6. Pull patella strap across knee and attach.

Print these instructions

669 East Industrial Park Drive - Manchester, NH 03109
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling.

If you have any questions about treatment for a specific condition, please check with a health care professional.