


45/100/125 kg



S.A.C.H. FOOT

- **S.A.C.H.** foot (Solid Ankle, Cushion Heel) is used for patients with low activity level (1 and 2).
- **S.A.C.H.** foot provides shock absorption at heel strike thanks to heel inserts of various hardnesses
- Step completion is ensured by the flexible forefoot.

SELECTION OF YOUR S.A.C.H. FOOT

Foot Shell		MEN				WOMEN & MEN		WOMEN		CHILDREN	
		With Toes				Without Toes		With Toes		With Toes	
Heel Height		10 mm		18 mm		18 mm		35 mm		12 mm	
Complete Item No for Order	Item N°	1A111		1A118		1A128		1A21		1A51	
	Side	-G	-D	-G	-D	-G	-D	-G	-D	-G	-D
	Foot Size (cm)	22	22	22	22	22	22	21	21	12	12
		23	23	23	23	23	23	22	22	14	14
		24	24	24	24	24	24	23	23	16	16
		25	25	25	25	25	25	24	24	17	17
		26	26	26	26	26	26	25	25	18	18
		27	27	27	27	27	27	26	26	20	20
		28	28	28	28	28	28				
		29	29	29	29	29	29				

S.A.C.H. feet come with a medium heel insert.

To optimize shock absorption at heel strike, the insert can be replaced by a softer or a harder one.

1A111 / 1A118 / 1A128			
Feet / Inserts	22 to 23 cm	24 to 27 cm	28 to 29 cm
Soft	1A11120-1S	1A11120-2S	1A11120-3S
Medium	1A11120-1N	1A11120-2N	1A11120-3N
Hard	1A11120-1F	1A11120-2F	1A11120-3F

1A21		
Feet / Inserts	21 to 24 cm	25 to 26 cm
Soft	1A11120-1S	1A11120-2S
Medium	1A11120-1N	1A11120-2N
Hard	1A11120-1F	1A11120-2F

UPPER LINK



• **1D111** p30
The combination of a **S.A.C.H.** foot and a **1D111** Multiaxis Ankle allows the patient to perform motions of plantar and dorsiflexion, as well as inversion-eversion.



• **1D43 / 1D44 / 1D48** p28
Pyramid Adapter for **S.A.C.H.** foot providing adjustment of equinus and valgus-varus.



• **1D38** p28
Pyramid Adapter for **1A51** child foot providing adjustment of equinus and valgus-varus.



• **1D05** p29
Fixed Tubular base.



• **1D07** p29
Tubular base providing adjustment of equinus.

• **1D60** (feet from 22 to 25) / **1D61** : (feet from 26 to 29) p29
Wooden base for exoskeletal prosthesis.

• **1D64** : (for **1A51** foot) p29
Wooden base for exoskeletal prosthesis.

RELATED PRODUCTS

• **1D20** p115
Junction plate ensuring a temporary attachment of the foot shell.



1A111 - Aesthetic Foot for Men, Sandal-Toe
Heel height : 10 mm

Item N°	Side	Size (cm)	Max. Patient's Weight
1A111	-G / -D	22 to 25	100 kg
1A111	-G / -D	26 to 29	125 kg

Approximate weight for Size 25 : 390 g



1A118 - Aesthetic Foot for Men, Sandal-Toe
Heel height : 18 mm

Item N°	Side	Size (cm)	Max. Patient's Weight
1A118	-G / -D	22 to 29	100 kg

Approximate weight for Size 25 : 410 g



1A128 - Foot for Men and Women, without toes
Heel height : 18 mm

Item N°	Side	Size (cm)	Max. Patient's Weight
1A128	-G / -D	22 to 29	100 kg

Approximate weight for Size 25 : 455 g



1A21 - Aesthetic Foot for Women
Heel height : 35 mm

Item N°	Side	Size (cm)	Max. Patient's Weight
1A21	-G / -D	21 to 26	100 kg

Approximate weight for Size 22 : 415 g



1A51 - Aesthetic Foot for Children
Heel height : 12 mm

Item N°	Side	Size (cm)	Max. Patient's Weight
1A51	-G / -D	12 / 14 / 16 17 / 18 / 20	45 kg

Approximate weight for Size 14 : 130 g