

- Pyramid Receiver for Bonding

This adapter will be bonded with epoxy (supplied) on a 30 mm Ø tube.



Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1D47	100 kg	Aluminium	85 g	52 mm	18 mm

- Threaded Pyramid insert M36 x 1,5 mm



Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K182	125 kg	Stainless steel	65 g	28 mm	6 mm

- Threaded Receiver insert M36 x 1,5 mm

Use with 1K179 or 1S145.



Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K183	150 kg	Stainless steel	75 g	23 mm	14 mm

- Threaded Pyramid

For screwing on
1P200-KD, 1P320-KD, 1P321-KD, 1P360-KD, 1P130-KD, 1P180-KD.



Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K184	125 kg	Stainless steel	94 g	30 mm	15 mm