

DOUBLE ADAPTER

Conversion table	
Kg	Lbs
45	99
80	176
100	220
125	275
150	330



• Double Pyramid Receiver

Avoid using a tube in tight spaces.

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K207-032	100 kg	Aluminium	85 g	32 mm	32 mm
1K207-045	100 kg	Aluminium	100 g	45 mm	45 mm
1K207-060	100 kg	Aluminium	115 g	60 mm	60 mm
1K207-075	100 kg	Aluminium	130 g	75 mm	75 mm



• Double Pyramid Receiver

Avoid using a tube in tight spaces.

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K207-HD-032	150 kg	Aluminium	100 g	32 mm	32 mm
1K207-HD-045	150 kg	Aluminium	130 g	45 mm	45 mm
1K207-HD-060	150 kg	Aluminium	145 g	60 mm	60 mm
1K207-HD-075	150 kg	Aluminium	160 g	75 mm	75 mm

• Telescopic Double Pyramid Adapter

Double adapter including two pyramid receivers. Can be adjusted in height by shortening, which allows the CPO to accurately adjust the prosthesis height. It can also provide rotation adjustment.

1K15001 and 1K155 : Threaded tube alone (in aluminium or steel)

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K150 (1)	100 kg	Aluminium	109 to 170 g	40 to 100 mm	40 to 100 mm
1K15001 (2)	125 kg	Aluminium	119 g	40 to 100 mm	40 to 100 mm
1K155 (2)	125 kg	Steel	119 g	40 to 100 mm	40 to 100 mm



Double adapter including a pyramid receiver and a pyramid. It allows for individual height adjustment of the prosthesis as its height can be shortened.

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K151	100 kg	Aluminium	141 to 201 g	44 to 89 mm	29 to 89 mm

5C413 Cutting Tool : Recommended to cut the 1K151 and 1K150 Adapters to the desired height.

• Double Pyramid Adapter

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K190	150 kg	Stainless steel	120 g	36 mm	7 mm

