

# NEW ICEROSS DERMO® UNIFORM LINERS

Iceross Dermo Uniform Locking



Iceross Dermo Uniform Cushion



## Changing from gel to silicone is easy with Iceross Dermo Uniform Liners.

Because now they are designed to easily retrofit the Uniform profiles commonly found in the market today, while still offering the popular features of our Dermo liners. Recommended for low to moderately active transtibial users, Iceross Dermo Uniform Liners are ideal for vascular patients and those with sensitive skin.

### Uniform Profile

- Retrofits common profiles without making a new socket
- 6 sizes accommodate most residual limbs

### Iceross Cover

- Elastic and smooth feeling material
- Easy donning/doffing
- Durable

### Seamless

- Seams cannot be felt through the silicone cushioning
- Welded distal ends on cushion liners

#### IMPACT LEVEL



**SUGGESTED BILLING CODES\***  
**L-5673, L-5679**

#### PRODUCT FEATURES



##### FABRIC COVER

The outer cover provides strength, durability, and protection against extreme wear and tear while allowing radial stretch and comfortable elasticity.



##### ACTIVE SKIN CARE

Vaseline® and Aloe Vera are incorporated into Dermosil® silicone to soothe, and are intended to maintain healthy skin.



##### SILKEN INNER SURFACE

Microengineered surface for superior skin conformance.




##### MATRIX

Unique integrated matrix allows radial stretch and limits vertical stretch for increased stability and skin comfort.



Life Without Limitations®

# NEW ICEROSS DERMO® UNIFORM LINERS



**Attention:** the current sizing has changed from the previous version. Please be sure to fit based on the above chart.

| 6mm      | Size | Proximal Circumference | Distal Circumference |
|----------|------|------------------------|----------------------|
| I-4X0615 | S    | 20-27cm                | 15-20cm              |
| I-4X0618 | M    | 25-37cm                | 18-23cm              |
| I-4X0620 | M+   | 36-52cm                | 20-26cm              |
| I-4X0625 | L    | 38-55cm                | 25-33cm              |
| I-4X0628 | L+   | 42-61cm                | 28-36cm              |
| I-4X0633 | XL   | 45-65cm                | 33-43cm              |

X = 1 (Locking) or 4 (Cushion)

The shape of the Dermo liners allows them to be used for retrofitting Alpha Locking and Cushion Liners.  
6mm liners can be used to retrofit 3mm liners in case of volume changes.

TEL (800) 233-6263  
FAX (800) 831-3160  
WEB OSSUR.COM

FOLLOW ÖSSUR ON



Life Without Limitations®