

Description	Measurement	Strap Lengths	Qty	Non-Sterile	Sterile
Arm/Leg	10" x 5"	--	1 Each	T5015NS	T5015ST
Hip/Lumbar	10½" x 11½"	--	1 Each	T5025NS	T5025ST
Knee/Shoulder	10½" x 11"	24" Detached (2) 36" Detached (1)	1 Each	T5035NS	T5035ST
		24" Detached (2) 36" Detached (1)	10/Case	T5035NS-10	T5035ST-10
Knee/Shoulder - Large	12½" x 11"	24" Attached (2)	1 Each	T5035NS-L	T5035ST-L
		24" Attached (2)	10/Case	T5035NS-L-10	T5035ST-L-10
Knee/Shoulder (w/ Attached Straps)	10½" x 11"	24" Attached (2)	1 Each	T5036NS	T5036ST
		24" Attached (2)	10/Case	T5036NS-10	T5036ST-10
Knee/Shoulder - Large (w/ Attached Straps)	12½" x 11"	24" Attached (2)	1 Each	T5036NS-L	T5036ST-L
		24" Attached (2)	10/Case	T5036NS-L-10	T5036ST-L-10
Shoulder/Knee (w/ Attached Straps)	10½" x 11"	36" Attached (1) 12" Attached (1)	1 Each	T5037NS	T5037ST
		36" Attached (1) 12" Attached (1)	10/Case	T5037NS-10	T5037ST-10
Shoulder/Knee - Large (w/ Attached Straps)	12½" x 11"	36" Attached (1) 12" Attached (1)	1 Each	T5037NS-L	T5037ST-L
		36" Attached (1) 12" Attached (1)	10/Case	T5037NS-L-10	T5037ST-L-10
Hand/Wrist	3½" x 5"	--	2 Each	T5055NS	T5055ST
Cervical	10" x 21"	--	1 Each	T5075NS	T5075ST
Foot/Ankle	10" x 11½"	--	1 Each	T5085NS	T5085ST