

T34 Hyperextension Orthosis

- Anatomically designed pubic pad and pelvic bar avoids pressure on pubis and prevents migration
- Sternal pad swivels to distribute pressure
- Latch-lock clamp closure allows easy donning and doffing
- Adjustable height and width
- Removable, washable pads

Indications: Osteoporosis, Thoracolumbar fractures T6-L1 and post surgical support

Contraindications:

- Three-column spinal fractures involving anterior, middle and posterior spinal structures.
- Compression fractures above T6, because segmental motion increases above the sternal pad.
- Ineffective in limiting lateral bending and rotation of the upper lumbar spine.



SPECIFICATIONS

Model No.	Size	Hip Circumference	Product Height
T34-02	Small	76-83 cm / 30-32½ in	43-47 cm / 17-18½ in
T34-03	Medium	83-94 cm / 32½-37 in	45-50 cm / 17½-19½ in
T34-04	Large	94-104 cm / 37-41 in	48-54 cm / 19-21 in
T34-05	X Large	104-114 cm / 41-45 in	53-59 cm / 21-23 in

Note: Original dimensions are in CM and conversions are rounded to the nearest ½".

Warranty: 90 days

USA
Tel: 1.800.492.1088
Email: o&porders@trulife.com

Canada
Tel: 1.800.267.2812
Email: customerservice@trulife.com

UK
Tel: +44 (0) 1142 618100
Email: info@trulife.co.uk

International/Ireland
Tel: +353 1 4511755
Email: info@trulife.com